CEMIX™: CEMENT MIXING SYSTEM

EN Instruction for use.....p.02

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
Ω Symbols.....p.35

CE 0123

Intended Use: device for mixing bone cement to be used in the vertebroplasty or kyphoplasty procedure

This device has to be used only in combination with medical devices intended for vertebroplasty procedure such as bone cement and syringe or injection devices EC certified and intended for vertebroplasty.

CHARACTERISTICS AND PERFORMANCE:

This device is a cement mixing system to manually mix the two components of the bone cement PMMA (powder and liquid components) during vertebroplasty or kyphoplasty procedure. The device is composed by a plastic container and a top part with a universal Luer-lock connection (for the connection with the syringes used in the procedure). The mixing of the two components (liquid and powder) of the bone cement is performed by shaking the device and by the action of the internal mixing ball

INDICATIONS:

The device is part of the Kyphoplasty systems manufactured by Biopsybell used to treat pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include haemangioma and giant cell tumour.

CONTRA - INDICATIONS:

The contra-indications of the kyphoplasty/vertebroplasty procedure are the followings:

- Infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Unstable vertebral fractures due to posterior involvement
- Previous damage to the pedicle wall (transpedicular access)
- Lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Retropulsing vertebral fragments with myelopathy
- A satisfactory response to conservative treatment
- Asymptomatic stable vertebral fractures

TARGET PATIENT:

The device is intended for the following patient population:

SEX: Male – Female

MATURATION OF THE SKELETAL SYSTEM: Skeletally mature patients

EXCLUSION: Newborns (<6 months) - Babies (>6 months to 2 years) - Skeletally immature patients

The device is used to treat patients with pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include haemangioma and giant cell tumour.

The patients to be excluded are the followings:

- Patients with infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Patients with vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Patients with unstable vertebral fractures due to posterior involvement
- Patients with previous damage to the pedicle wall (transpedicular access)
- Patients with lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Patients with retropulsing vertebral fragments with myelopathy
- Patients with a satisfactory response to conservative treatment
- Patients with asymptomatic stable vertebral fractures

INTENDED USERS AND ENVIRONMENT USE CONDITIONS:

The device is for professional use and, in particular, must be used by qualified and trained radiologists, neuro-radiologists, interventional radiologists, orthopaedic surgeons and neuro surgeons, pain therapists. The professional must have documented basic skills about kyphoplasty. The device is to be used inside a sterile surgical field such as operating and angiographic rooms.

WARNINGS, PRECAUTIONS AND SIDE-EFFECTS:

- Do not use the instrument if the package is damaged.
- Do not use the instrument beyond the expiry date.
- The device is single-use and must be destroyed after use. Reuse is strictly prohibited. Reuse or resterilization may compromise the structural integrity of the device and/or cause damage to the device which may cause injury, disease or death of the patient. Reuse or resterilization may also create a risk of device contamination and/or cause infection of the patient or cross-infection, including but not limited to transmission of infectious diseases from one patient to the next. Device contamination may cause harm, disease or death of the patient.
- Do not reuse a component more than once.
- The system may exclusively be used by physicians. After use, the device may constitute a risk of contamination and/or infection. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.
- In order to minimize the risks or other complications associated with use of the cement for vertebroplasty/kyphoplasty, the following is recommended:
- For all patients that have acute osteoporotic fractures of the vertebral body, a period of conventional treatment should be considered.
- Only qualified and appropriately trained physicians should perform the vertebroplasty/kyphoplasty procedure.
- Use only cements intended for vertebroplasty/kyphoplasty and carefully read the instructions for use of the cement.
- It is strictly prohibited to modify, restructure, repair or resterilise the device for reuse.



INSTRUCTION FOR USE

Preliminary operations:


Before use check that the package is not damaged and that the set and its components are not damaged or kinked; if so, inform the manufacturer or your area representative. Take the devices out of the package operating in aseptic conditions.

INSTRUCTIONS:


1. Unscrew the Top part of the mixer.
2. Close the system with Top part provided (screw it with the container) and shake manually very vigorously

3. Place inside the mixer the Powder and then the Liquid that will form the cement for vertebroplasty (this operation must be performed as quickly as possible).

4. After making sure the cement has the right consistency (following the instructions for use of the bone cement manufacturer), unscrew the Luer Lock Cap and connect to the side a specific syringe with Luer- Lock connection (or a specific injection device) and proceed to suck up the cement according to the instructions for use.

5. Once the necessary quantity has been sucked up, disconnect the syringe from the mixer.
6. After terminating the use dispose of the device according to applicable laws and medical practice.

PACKAGING, STERILIZATION, SINGLE USE:


The system is individually packed in double pouch made of medical grade paper 60g. and transparent film PET/PP. Do not use the device if the package is open or damaged or contaminated. It is sold in packs of 5 pieces. The needle is sterilized by Ethylene Oxide (EtO) in compliance with ISO 11135 requirements. The device is not reusable. The device is single-use and must be destroyed after use.

DISPOSAL:

After use, the device may constitute a risk of contamination and/or infection and puncture or injuries caused by sharp parts of the device. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.

NOTICE FOR USER AND/OR PATIENT: ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE SHOULD BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED.

Biopsybell declines all responsibility for any damage caused by improper use or different from what is indicated in this instruction leaflet.

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