
EASYNJECT™

EN

Bone cement delivery system.....p. 02

IT

Sistema di iniezione del cemento osseo.....p. 20

BG

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
Intended Use: device for injection of bone cement during vertebroplasty or kyphoplasty procedure.

CHARACTERISTICS AND PERFORMANCE:

The device is composed by a syringe cylinder in which is contained the bone cement to be injected and a threaded plunger for pushing out the cement in the injection procedure. The device is equipped with a rear knob that allows to retract and advance the plunger. The device has a rotating luer lock connection for connecting the extension tube or a mixing device used in the procedure.

Compatible with low-, medium- and high-viscosity cements.

- Ergonomic grip.
- 10ml single-use syringe included.
- Millimetric plunger advance that allows accurate and controlled cement injection.
- Allows fast and direct injection by pressing the button on the grip and moving the plunger forward.
- Equipped with a 30cm extension tube.

This device has to be used only in combination with medical devices intended for vertebroplasty procedure such as bone cement EC certified and intended for vertebroplasty.

INDICATIONS:

The device is part of the Kyphoplasty systems manufactured by Biopsybell used to treat pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include hemangioma and giant cell tumor.

CONTRA-INDICATIONS:

The contra-indications of the kyphoplasty/vertebroplasty procedure are the followings:

- Infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Unstable vertebral fractures due to posterior involvement
- Previous damage to the pedicle wall (transpedicular access)
- Lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Retropulsing vertebral fragments with myelopathy
- A satisfactory response to conservative treatment
- Asymptomatic stable vertebral fractures

Target Patient:

The device is intended for the following patient population:

SEX: Male – Female

MATURATION OF THE SKELETAL SYSTEM: Skeletally mature patients

EXCLUSION: Newborns (<6 months) - Babies (>6 months to 2 years) - Skeletally immature patients
The device is used to treat patients with pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include hemangioma and giant cell tumor.

The patients to be excluded are the followings:

- Patients with infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Patients with vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Patients with unstable vertebral fractures due to posterior involvement
- Patients with previous damage to the pedicle wall (transpedicular access)
- Patients with lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Patients with retropulsing vertebral fragments with myelopathy
- Patients with a satisfactory response to conservative treatment
- Patients with asymptomatic stable vertebral fractures

INTENDED USERS AND ENVIRONMENT USE CONDITIONS:

The device is for professional use and, in particular, must be used by qualified and trained radiologists, neuro-radiologists, interventional radiologists, orthopaedic surgeons and neuro surgeons, pain therapists.

The professional must have documented basic skills about kyphoplasty.

The device is to be used inside a sterile surgical field such as operating and angiographic rooms.

WARNINGS, PRECAUTIONS AND SIDE EFFECTS:

- Do not use the instrument if the package is damaged.
- Do not use the instrument beyond the expiry date.
- The device is single-use and must be destroyed after use. Reuse is strictly prohibited. Reuse or resterilization may compromise the structural integrity of the device and/or cause damage to the device which may cause injury, disease or death of the patient. Reuse or resterilization may also create a risk of device contamination and/or cause infection of the patient or cross-infection, including but not limited to transmission of infectious diseases from one patient to the next. Device contamination may cause harm, disease or death of the patient.
- Do not reuse a component more than once.
- The system may exclusively be used by physicians. After use, the device may constitute a risk of contamination and/or infection. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.
- In order to minimise the risks or other complications associated with use of the cement for vertebroplasty/kyphoplasty, the following is recommended:
- For all patients that have acute osteoporotic fractures of the vertebral body, a period of conventional treatment should be considered.
- Only qualified and appropriately trained physicians should perform the vertebroplasty/kyphoplasty procedure.
- Use only cements intended for vertebroplasty/kyphoplasty and carefully read the instructions for use of the cement.
- Monitor the procedure using high-quality imaging devices in order to detect PMMA leakage or seepage.
- Closely monitor the patient's blood pressure during and immediately after the procedure; multiple treatments may increase the risk of sudden blood pressure drops due to release of PMMA monomer into circulation. Not more than 3 vertebra should simultaneously be treated in one session.
- Careful diagnosis should be performed and particular precautions taken if the procedure is performed on patients with spinal tumours that have eroded the posterior wall of the vertebral body.
- It is strictly prohibited to modify, restructure, repair or resterilize the device for reuse.

INSTRUCTION FOR USE:

Preliminary operations

Before use check that the package is not damaged and that the set and its components are not damaged or kinked; if so, inform the manufacturer or your area representative. Take the devices out of the package operating in aseptic conditions.

Instructions

1. Prepare and mix the bone cement used for the procedure according to the manufacturer's instructions.
2. Make sure that the safety lock - active when the plunger is fully inserted and screwed in up to the end stop - is released. To release the lock, turn the Rear Knob anticlockwise until the red O-ring is positioned approximately on the 1ml notch of the syringe.
3. Connect the rotating luer lock to the extension line or to the cement mixing device used to prepare the bone cement.
4. Hold down the Release button located on the side of the grip (release button) and at the same time pull the threaded plunger backward to allow cement aspiration into the syringe cylinder.
5. Gradually move the threaded plunger forward in the syringe by turning the rear knob clockwise in order to inject the cement.
6. To stop cement injection, turn the rear knob anticlockwise by 180°.

PACKAGING, STERILIZATION, SINGLE USE:

The system is individually packed in double pouch made of medical grade paper 60g. and transparent film PET/PP. Do not use the device if the package is open or damaged or contaminated. It is sold in packs of 2 pieces. The device is sterilized by Ethylene Oxide (EtO) in compliance with ISO 11135 requirements. The device is not reusable. The device is single-use and must be destroyed after use.

DISPOSAL:

After use, the device may constitute a risk of contamination and/or infection and puncture or injuries caused by sharp parts of the device. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.

NOTICE FOR USER AND/OR PATIENT: ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE SHOULD BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED.

Biopsybell declines all responsibility for any damage caused by improper use or different from what is indicated in this instruction leaflet.

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